



KF

Karl Fischer Coulometer



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Karl Fischer Coulometer

DSKFT-20V



GENERAL FEATURES

- LCD display screen.
- Generator cell with diaphragm.
- Selectable units including μg , mg, %, mg/L, g/L, etc.
- Settable parameters, including measurement unit, polarization current, generator current, endpoint potential, stop criterion and etc.
- Data can be stored and easily transferred to printer via RS-232 communication interface.
- Reset feature automatically resumes all settings back to factory default options.

INCLUDES

- JB-11 stirrer
- Electrolysis electrode
- Work electrode



Model	DSKFT-20V	
Coulometric Ttration	Range	10 μg ~20mg
	Polarization Current Accuracy	$\pm 0.2\mu\text{A}$
	Working Current	10mA, 20 mA, 50 mA, 100 mA
	Working Current Accuracy	$\pm 0.5\%$
	Accuracy	$\pm (5\%+3) \mu\text{g}$
	Repeatability	$\leq 3\%$
General	Power Supply	AC Adapter, 200-240 VAC input
	Dimensions	300 $\times$ 235 $\times$ 100 mm
	Net Weight	2 kg
	Dimensions (Carton)	370 $\times$ 305 $\times$ 330 mm
	Shipping Weight	3 kg