



MP Microplate Centrifuge



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

FEATURES



**Rotor for microplate
more working
capacity**



**Brushless motor
with strong power**



**Speed adjustable
and timing**



LCD display

Darsun Scientific

Microplate Centrifuge

DSL-500-MP



Model	DSL-500-MP
Max Capacity	96-well PCR microplate * 4
Speed Range	500 - 2000 rpm
Max RCF	500xg
Semidiameter Of Rotor	10.2cm
Protection Function	Lid opening protection, cavity over temperature protection, vibration projection
Timing	30s - 59min50s
Power	50W
Power Supply	AC 110V - 220V 50/60Hz