



MC

Mini Centrifuge



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific

FEATURES



Stop running after opening lid



Low noise

Darsun Scientific

Mini
Centrifuge

DSMC-4K/DSMC-7K/DSMC-10K/
DSMC-12K



Model	DSMC-4K	DSMC-7K	DSMC-10K	DSMC-12K
Max Speed	4000rpm	7000rpm	10000rpm	12000rpm
Max RCF	980xg	3010xg	6149xg	8855xg
Rotor Capacity	0.2ml*16 + 1.5 / 2.0ml*8		1.5 / 2.0ml*8	
Power	15W	20W	30W	40W
Dimensions (mm)	160*135*110	160*135*110	160*135*110	160*135*110
Weight (kg)	0.8	0.8	0.8	0.8