



MM

Multiple Magnetic Stirrer With Heating



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Multiple Magnetic Stirrer With Heating

DSMHS4P/DSMHS6P



FEATURES



Heat more uniformly



LCD display



304 stainless steel plate,
Corrosion resistant



Durable brushless motor

Model	DSMHS4P	DSMHS6P
Control	Asynchronous control	
Number Of Stirring Points	4	6
Max Capacity	4*500ml	6*500ml
Working Plate Size	4-(110*110mm)	6-(110*110mm)
RPM	50-1200 rpm	
Plate Material	304 Stainless steel	
Plate Temperature	120°C	
Timing	1min-19h59min	
Bar Size (mm)	Z30	Z30
Power	240W	360W
Voltage	AC 220V	