



OS Orbital Shaker



Darsun Scientific Inc, USA
503 N Gould St, Ste R, Sheridan, Wyoming
82801, United States

Email : info@darsunscientific.com
Phone : +1-3078722810
Web : www.darsunscientific.com

DarsunScientific

Darsun Scientific

Orbital Shaker

DSOS-40



Product Features

- 1.Speed: 40-200rpm
- 2.Digital panel with speed and time control
- 3.Working pan size 32X32cm,with max load 5KG
- 4.Continuous operation and timed operation range 1min-23h59min
- 5.DC brushless motor, safe, long life and free of maintenance. Operating silent and stable

MODEL NO	DSOS-40
Motor type	Brushless motor
Plate size	320*320
Turn	Circular
Maximum load	5 kg
Orbit diameter	20mm
Speed range	40-200 rpm
Timing range	1min-23hr59min
Voltage	100-220V,50/60HZ
Allowable	5-40°C 80% RH
Size	435*330*140mm
Net weight	8.1 kg