

DEEP FREEZER

300/400L, -10~-25°C



DEEP FREEZER 300/400L, -10~-25°C

Efficient compressor, accurate temperature control

Suitable for hospital, blood-bank, center for Disease Prevention, scientific research institution, animal husbandry industry, electronic lab and chemical lab, etc

Product Features

TEMPERATURE CONTROL SYSTEM

Controlled by micro-computer, digital display of inner temperature, the inner temperature can be adjusted at the range of -18°C~-25°C; Temperature control system, you can set the temperature alarm according to requirements.

SECURITY CONTROL SYSTEM

Three malfunction alarms: high temperature alarm, low temperature alarm, and sensor malfunction alarm;
Two alarming methods: buzzer alarming and digital flashing alarming;
Protection function: power-on delay, time between stops

REFRIGERATION SYSTEM

Famous brand compressor and fan motor;
Thickened temperature preservation layer, ultra-microcellular foaming technology, outstanding temperature preservation;
CFC-free foaming, CFC-free refrigerant, more environmental friendly;
Optimized evaporation system and condensation system, strong cooling, low energy consumption

HUMANIZED DESIGN

Drawer design, convenient for deposit;
Wide voltage design;
Caster design at the bottom, free to move and easy to fix;
Wide climate zone design, suitable for the ambience with temperature of 10°C~32°C.
Safety door lock, preventing the door opened by chance.
Pressure-balanced valve design, easy to open the door

Model No.	Total Valid Volume (L)	Inner Temperature (°C)	Rated Voltage/Frequency	External Dimension (W×D×H) (mm)	Internal Dimension (W×D×H) (mm)
DS-25L300	300	-18~-25	~220V/50	600×625×1875	440×448×1673
DS-25L400	400	-18~-25	~220V/50	775×720×1900	605×529×1699

Cooling Curve

