



WB

Water Bath



Darsun Scientific Inc., USA
30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

Email: Info@darsunscientificinc.com

Phone: +1-307-427-2810

Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Water Bath

DS-1-WB100/DS-2-WB100
DS-4-WB100/DS-6-WB100
DS-8F-WB100



FEATURES



One-piece stainless steel tank



High corrosion resistance



High precise temperature control



Digital display thermostatic

Model	DS-1-WB100	DS-2-WB100	DS-4-WB100	DS-6-WB100	DS-8F-WB100
Temperature Range	Room temp-100℃				
Power	300W	600W	800W	1200W	1500W
Accuracy	±0.1℃				
Working Chamber Size (mm)	160*160*130	300*150*140	300*300*140	450*300*130	620*300*120