



WB

Water Bath With Magnetic Stirrer



Darsun Scientific Inc., USA
30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

Email: Info@darsunscientificinc.com

Phone: +1-307-427-2810

Web: www.darsunscientificinc.com



Darsun Scientific

Water Bath With Magnetic Stirrer

DS-WB-6S/DS-WBS-6P



FEATURES



Brushless motor



Strong Stirring Power



PID temperature control



Stainless steel inner tank

Model	DS-WB-6S	DS-WBS-6P
Speed / Increment	/	100-1200RPM / 50RPM
Temperature Range	Room Temp. +5°C -100°C	Room Temp. +5°C -100°C
Liner Size	200*200*130mm(depth)-2L flask or 3L beaker	
Product Size (mm)	260*265*325	260*265*325
Display Accuracy	0.1°C	0.1°C
Power	800W	
Timing Range	1min-23h59min	