



OS

3D Oscillator Shaker



Darsun Scientific Inc, USA
30 N Gould St, Ste R, Sheridan, Wyoming
82801, United States

Email: info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

3D Oscillator Shaker

DSSK-3D



FEATURES



Timing function



Brushless motor with strong power



3D mixing



Suitable for a variety of containers

Model	DSSK-3D
Type	3D Oscillating
Diaplay	LCD
Tit Angle	10°
Max Load	0.8kg
Speed	10 - 80rpm
Timing	1min - 8h
Allowed Ambient Humidity Temp.	5~40°C 80%RH
Voltage	AC 100-240V, 50/60Hz
Power	10W
Dimension (mm)	184*184*130