



Chemiluminescence immunoassay Analyzer



Darsun Scientific Inc, USA

30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific chemiluminescence immunoassay Analyzer



Product Features

Product Type: Bench-Top Fully automatic chemiluminescence immunoassay Analyzer

Product weight: Net weight: 62kg

Processor: RAM 2GB, main frequency 4x 1GHz or above

Security Interface: USB 2.0 & Ethernet RJ45 interface

1. Total 90 test parameters including infectious markers (HIV, HCV, HBV & Syphilis), Dengue range (NS1 Ag, IgM & IgG), cardiac markers, fertility markers, etc.
2. Solid Phase Carrier: Magnetic Particles.
3. Processing mode: Batch & continuous random access.
4. Result Throughput is 120 tests/hour (max.), first result within 15 min.
5. 30 sample position & 10 reagent positions to process multiple test samples at a time.
6. Patented Single steel S/R Probe (needle) for sample & reagent addition contributing in high reliability & more accuracy in test results.
7. Single assay cup for reagent & sample processing leading to lower cost in consumables.
9. Large data storage capacity for storing around 5,00,000 patient data. HIS & LIS compatible with application support.

DARSUN
SCIENTIFIC