



AH

Automatic Hematology Analyzer



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Automatic Hematology Analyzer

DSK5000P



DARSUN
SCIENTIFIC

- 3-part differentiation of WBC •Double Chambers for WBC and RBC/PLT
- Throughput: 60 samples per hour •Intelligent reagent management function: Reagent information and residual display •Strong data management system •Multilingual for option •LCD with high resolution
- Bar code scanner (optional) •Reagent alarm system (optional)

21Parameter + 3Histogram WBC, LYM#, MID#, GRA#, LY%, MTD%, GR%, RBC, HGB,HCT, MCV, MCH, MCHC, RDW-SD, RDW-CV, PLT, MPV, PDW, PCT, P-LCR,P-LCCandHistogramforWBC,RBC,PLT

Model	DSK5000P
Predilute	20 μ L
Whole blood	10 μ L
Printout	Built-in thermal printer
Temperature	10-35 $^{\circ}$ C
Humidity	$\leq$ 90%
Display	8" color touch screen display
Menu	Review, Test, Maintenance, QC, Setup, Shutdown, Prime reagents, System diagnosis, Aperture rinse, Soak and Rinse, Calibration, instrument config
Throughput	60 samples per hour
Input / Output	1 RS232, 1 keyboard , Mouse; 4 USB, LAN
Power Requirement	Power Requirement
Aperture Diameter	WBC 80 μ M RBC/PLT 70 μ M Carryover WBC $\leq$ 1.5%,RBC,HGB
Dimension	Instrument:- 440mm(H) x 350mm(W) x 440mm(D) Net:17kg Package:-570mm(H) x 520mm(W) x 600mm(D) Gross:21.8kg

DarsunScientific