



Automatic coagulometer



Darsun Scientific Inc, USA
30 N Gould St, Ste R Sheridan, Wyoming
82801 United States

Email: info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific Automatic coagulometer

DS-202



Test items Main test:

FIB Others:

Hep Proc Pros

Coagulation factor II,V,VII,VIII,IX,X,XI,XII

AT-III

AT Rep T RVVT RVVCT RVVRT LA

Product Features

Double magnetic route, stable performance , not liable to interference. 2. 7" Color touch screen 3. 100 assay parameters and 1000 test results can be saved 5. 58mm thermal printer 7. Clock CMOS chip inserted 8. Standard 232 serial interface 9. Power supply: 220V±15%, 50Hz, 80VA

Test method: double magnetic route 4. Temperature: 37°C±0.5°C 5. adding sample: by a timing linkage 6. precision: CV ≤3% 7. speed: ≤300s(depend on setting)

Model	DS-202
Operation	Automatic
Number of channels	4-channel
Test type	PT, APTT, TT
Configuration	Benchtop
Other characteristics	With touchscreen
Sample volume	0.2 ml (0.006763 US fl oz)
Number of sample positions	8-unit, 16 unit
Reagent volume	Min.: 0.05 ml (0 US fl oz) Max.: 0.2 ml (0.01 US fl oz)
Number of reagent positions	6 unit
Weight	5 kg (11 lb)
Length	370 mm (14.6 in)
Width	290 mm (11.4 in)
Height	120 mm (4.7 in)