



## Automatic veterinary 3-diff Hematology Analyzer



**Darsun Scientific Inc, USA**  
30N Gould St Ste R Sheridan, Wyoming  
82801, United States

**Email:** [Info@darsunscientificinc.com](mailto:Info@darsunscientificinc.com)  
**Phone:** +1-3074272810  
**Web:** [www.darsunscientificinc.com](http://www.darsunscientificinc.com)

**Darsun**Scientific



## DarsunScientific

### Automatic veterinary 3-diff Hematology Analyzer

#### DSK-VET3



| Model                | DSK-VET3                                                                                              |
|----------------------|-------------------------------------------------------------------------------------------------------|
| 21 parameters        | WBC LYM#,MID#, GR# LYM% MID% ,GR% ,RBC,HGB,HCT,MCV,MCH RDW-CV, RDW-SD, PLT,MPV, PDW, PCT, P-LCC,P-LCR |
| 3 Histograms         | WBC,RBC,PLT                                                                                           |
| Species              | 12 preset species including Canine ,Equine , Bovine, Murine , Swine , Lipoid 4 optical species        |
| Sample volume        | Whole blood 10 µl                                                                                     |
| Storage capacity     | Auto-save , more than 200.000 analysis results available                                              |
| Throughout           | 60 tests/hour                                                                                         |
| Fate search mode     | By sample ID, patient ID, date, multi-keywords. Loading/extracting data, offline reading              |
| Display              | 10.4 inch colorful touch screen                                                                       |
| Operation mode       | L-J ,X-R,X-B, etc.                                                                                    |
| Printout             | Thermal print (defaulted ) available for external printers optional printing templates and format     |
| Data communication   | Wifi, Bluetooth                                                                                       |
| Ports                | LAN,X,RS232X1,HDMIX1,USB port x4                                                                      |
| Language             | English                                                                                               |
| Optional Accessories | Handheld barcode scanner                                                                              |
| Dimension & weights  | 400x350x440mm 17 kg                                                                                   |
| Power supply         | 110-240V , 50/60Hz                                                                                    |

110-240V , 50/60Hz