



Basic Power Supply



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Basic Power Supply

DS-300



Product Overview

DS series Basic power supply comprises 2 models, including :
DS-300 and DS-600 , which can meet the conventional horizontal electrophoresis, small and medium-sized vertical electrophoresis, membrane acetate electrophoresis experiment requirement , the lightweight and compact, simple and practical characteristics, make it popular in for basic teaching.

Product Features

Type of Output: Constant voltage or constant current Output range: DS-300: 5-300V, 1-300mA, maximum 90W Increment: 1V, 1mA
Timing Range: 1min ~ 9hr, 59min Display: LED
Output jacks: Two sets of output jacks Dimension (LxWxH): 280x 237x 118 (mm) Weight: 3.0 (kg)
Molding shell, touch keys, microprocessor intelligent control;
Automatic memory;
Standard timing operation;
Automatic detection of no-load, over-load

