



Blood Coagulation Analyzer



Darsun Scientific Inc, USA

30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Blood Coagulation Analyzer

DS-CA02/DS-CA04



The blood coagulation analyzer adopts the principle of optical colorimetry for detection. After the reagent is mixed with the plasma, fibrinogen is converted into fibrin and coagulated, which leads to the change of the optical density of the test specimen, and the instrument detects the coagulation end point.

1. Dual-channel independent detection, can analyze different items at the same time
5-inch high-definition touch screen, easy to operate
2. Open system, good compatibility, supports reagents of various specifications
200,000 results storage
3. Equipped with multiple pre-temperature positions to improve detection efficiency
Built-in thermal printer
4. The pipette automatically triggers the measurement, and the measurement time can be accurate to tenths of a second
5. Less reagent consumption

Model	DS-CA02	DS-CA04
Test Principle	Optical colorimetry	
Test Items	PT, APTT, TT, FIB	
Test Channel	2 (can test different projects at the same time)	4 (can test different projects at the same time)
Wavelength	470nm	
Detection Time	Normal sample detection time is 20 to 40 seconds, can be set arbitrarily from 20 to 99 seconds	
Sample Incubation Position	12	24
Reagent Incubation Position	5	6
Detection Position Temperature	37.0 ±1.0°C	
Incubation Position Temperature	37.0 ±1.0°C	
Incubation Time	1 ~ 999s	
Reaction Time	Minimum dosage 20μl	
Reagent Dosage	20μl ~ 40μl	
Sample Volume	5-inch high-definition color touchscreen	
Display	20,000 results storage, auto-save on power off	
Storage	RS232 serial interface, supports LIS system	
Communication Interface	Built-in thermal printer	
Printer	AC100 ~ 240V, 50/60Hz	
Power Supply	70VA	
Consumption	10°C ~ 30°C	
Ambient Temperature	≤80%	
Relative Humidity	86kPa ~ 106kPa	
Atmospheric Pressure	360 ×210 ×120 mm	
External Size (W×D×H)	3.5 kg	
Net Weight	376 ×246 ×316 mm	
	4.96 kg	
	390 ×250 ×135 mm	
	445 ×305 ×265 mm	
	7 kg	