



Constant Temperature Circulator



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Constant Temperature Circulator DS-1050



ProductFeatures

DS-1050 Constant Temperature Circulator used for control keeping experiment's temperature at stable values during electrophoresis run s, DS-1050 also can be connected with the tank or other lab equipment to serve as a constant temperature controller.

Model	DS-1050
Temperature control range	-5~35 (°C)
Temperature accuracy	+ 0.1 (°C)
Circulating fluidrate (L/ min)	5 (L/MIN)
Tank volume (ml)	2.2 (L)
Power	300 (W)
Dimension (LxWxH)	196x406x482 (mm)
Weight	20.0 (kg)

Partsname	Quant it y
M ai nfram e	1
Powercord	1