




Constant temperature shaking incubator



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Constant temperature shaking incubator



Product Name	Constant temperature shaking incubator
Oscillation frequency	10–300 RPM (top layer 10–250 RPM)
Oscillation frequency accuracy	±1 rpm
Shaking plate amplitude	26mm (standard) 50mm (optional)
Temperature control range	4 to 60 °C
Temperature regulation accuracy	Plus or minus 0.1 °C
Temperature uniformity	±0.3°C (at 37°C)
Display mode	LCD (touchscreen)
Control methods	P.I.D Microcomputer intelligent control
Timing range	0–999 hours (multiple Settings available)
Shake plate size (L × W)	700mm * 430mm
Volume	200 liters
Standard configuration	Universal spring clamps (fixed clamps optional)
Power source	AC220 ± 10% 50 to 60Hz
External dimensions (L × W × H)	1 floor 1080mm × 685mm × 650mm (including base height)
Weight	Single layer 180kg
Opening method	Floating top opening
Power	1500W