



DB

Dry Bath Shaking Incubator



Darsun Scientific Inc, USA
30 N Gould St, Ste R, Sheridan, Wyoming
82801, United States

Email: info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Dry Bath Shaking Incubator

DS-DB100C/DS-DB100



FEATURES



LCD display



Over temperature protection



Program memory, 9 temperature points



Two modules, freely combined

Model	DS-DBS100C	DS-DBS100
Temperature Range	0°C ~ 100°C	15°C ~ 100°C
Temperature Stability	±0.5°C	±0.5°C
Display Accuracy	0.1°C	0.1°C
Heating Rate	5.5°C/min	5.5°C/min
Cooling Rate	2.5°C/min(100°C-Room temperature) 0.5°C/min(Below room temperature)	/
Speed / Timing	200-1500rpm / 1s-23h59min	
Mixing Diameter	3mm	
Program Memory	9 grades	
Power	200W	
Voltage	100-240V 50/60HZ	
Size	200*235*120mm	
module options	Numbers and aperture of test tube	
DBs1	0.2ml*96 96-hole PCR microplate	
DBs2	1.5ml *24 centrifuge tube	
DBs3	2ml *24 centrifuge tube	
DBs4	15ml *8 centrifuge tube	
DBs5	50ml *4 centrifuge tube	
oem of modules	Modules can be customized according to customer needs	
DBS-T	Low module support for thermal cover (thermal corver and transpatent house cannot be used at the same time)	