



Fluorometer



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Fluorometer

DS - FM201



Product Overview

The fluorometer is a small, portable fluorescence quantitative detector used for high-precision fluorescent nucleic acid quantitative determination of DNA, RNA, micro-RNA and protein concentrations. It has two channels and can quickly analyze two fluorescent signals of one sample.

Product Features

5-inch color touch screen with android system, friendly user interface. Quickly and accurately quantify DNA, RNA and proteins in 3 seconds. High sensitivity, the lowest detection limit can reach 0.5pg/ μ l with dual-color fluorescence channels for quantification of nucleic acids and proteins, measuring two different fluorescence in one analysis. Compatible with common reagents on the market. Can store more than 100,000 sample results and output data through USB port. High quantitative measurement accuracy can be achieved with only 1~20 μ l sample. Support multiple operating languages. 5 orders of magnitude linear dynamic range.

Model	DS-FM-201
Sample Throughput	1
Sample Tube	0.5 ml PCR tube (thin-walled tube)
Display	5-inch color touch screen
Data Output	USB/WIFI
Light Source	Blue 460~650nm
Excitation Filter	Blue 460~480nm Red 630~650nm
Emission Filter	Green 500~535nm Red 670~710nm
	1~20 μ l
Sample Volume	3s
Processing Speed	L200 x W100 x H54 mm
Dimension	0.43 kg
Weight	100~240V; 50/60Hz

Power	Supply	Common Fluorescent reagents	Application
470nm		Pico green oligreen RibogreenGfp,protein	Ds DNA ssDNA quantification green fluorescent protein GFP Detection prorenin quantification Rhodamine detection
525nm (optical)		fluorescein quant-it Rhodamine Cy3 RFP VybrantCytotoxicity	Cy-3 fluorescent labeling/red fluorescent protein RFP gatedetection RNA quantification Cy5 fluorescent labeling detection RNA quantification
650nm		Cy5 quant-it RNA	