



WB

Infrared Sterilizer



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Infrared Sterilizer

DS-900/DS-900S



FEATURES



One-piece stainless steel tank



High corrosion resistance



No open-flame, more safe



Easy to clean

Model	DS-900	DS-900S
Power	150W	240W
Max. Temperature Inside The Heating Hole	825°C±50°C	
Heating Temperature	Two temperatures available 400°C /850°C	
Max. Diameter Of Sterilizer Goods (mm)	Φ14	Φ35
Length Of Heating Zone (mm)	150	100
Warm-up Time	10min	