



LC

Large Capacity Magnetic Stirrer



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Large Capacity Magnetic Stirrer

DSMS-5L/DSMS-10L
DSMS-15L/DSMS-30L
DSMS-50L



FEATURES



Timing function



Knob operation



Forward and reverse



Large capacity

Model	DSMS-5L	DSMS-10L	DSMS-15L	DSMS-30L	DSMS-50L
Control Mode	Display & Knob				
Display Method	LCD display				
Power	10W	20W	20W	30W	50W
Working Plate (mm)	180*180	220*220	220*220	260*260	330*330
Max Stirring Volume	5L	10L	15L	30L	50L
Speed Range	50-1500rpm	50-1500rpm	50-1500rpm	50-1300rpm	50-1100rpm
Speed Regulation Accuracy	5rpm				
Timing Range	1min-23h59min				
Stir Bar Size (mm)	40	50	60	80	100
Product Size (mm)	280*180*65	315*220*65	315*220*65	350*260*85	435*330*85