



Manual colony counter



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Manual colony counter

DS-A11



Product Overview

- 1 Displays average values
- Bright LED lighting
- 2 Three-digit display for easy readings Adjustable pressure
- Beeping noise prevents errors
- 3
-
- 4

Product Features

Manual colony counter DS-A11 accommodates a wide range of colony counts from 0 to 999 for accurate measurements. This colony counter fits petri dishes from 50 to 150 mm, offering versatility. It features multi-mode functionality for simultaneous counting from two petri dishes. Equipped with a pressure-sensitive touch-tip pen marker. Ensures thorough examination and accurate counting of colonies through cumulative marking. Provides nine times magnifications for easier counting.

Model	DS-A11
Counting range	0 to 999
Compatible petridish	50 to 150 mm
Magnification	3 to 9 times
Digital display	3 digits
Rack	2
Relative Moisture	0.8
Ambient temperature	5 to 50°C
Protective class	IP21
Power consumption	20 W
Power	AC 100 to 240V/ 50 to 60 Hz
Dimension	360x300x180mm
Weight	4.0 kg