



Microplate Centrifuge



Darsun Scientific Inc, USA
30N Gould St Ste RSheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Microplate Centrifuge

DS-450



Product Overview

The DS-450 bench top micro plate centrifuge seamlessly integrates cutting-edge technology with substantial capacity.

Product Features

- LED display for easy parameter adjustment during operation, featuring a user-friendly program design. Customizable modes such as short spin, timed cycles, and rotational speeds ensure efficient and stable tests.
- Impulse centrifugation provides instant acceleration to the desired speed, with automatic braking and fast shutdown for added safety.
- Equipped with a maintenance-free DC brushless motor, allowing rapid acceleration and deceleration, enabling quick and efficient centrifugation of multiple samples.
- Compact and uniquely designed, the centrifuge is user-friendly, allowing easy handling of microplates. The lid opens up to 90° for convenient access.
- Efficiently disperses hanging liquids on tube walls before or after PCR tests, significantly enhancing test accuracy.
- Manufactured in compliance with international safety regulations,
- Safety features include an electrical lid interlock preventing door opening while the rotor is running. The machine will not operate if the lid is open, displaying an error code.
- Emergency lid-lock release function activates in case of a power failure.
- Equipped with a gas hinge for lid drop protection.
- Various safety protections, including over speed mechanisms, are integrated to ensure user and machine safety.

MODEL NO		DS-450			
Speed range(rpm)		0-3000rpm adjustable			
Max Capacity(ml)		2x96 well PCR plate / microplate, 24x8x0.2ml PCR strip tube, 192x0.2ml PCR tube			
Display		LCD			
Dimension(mm)		230x195x195mm			
Speed accuracy		±10rpm			
RCF range (xg)		0-600xg adjustable			
Time range		Short spin/ 1~99 seconds/ 1~99 minutes			
Motor		Frequency Converter motor			
Power supply		DC 24V power adapter (AC 100~230V 50/60Hz 5A), 55W			
Net weight(without rotor)		3.9 KGS			
No.	Item	Capacity	Rotational speed	Max.RCF	Note
NO.1	Adapter	24×8×0.2ml	3000r/min	600xg	PCR tube-strip or 0.2ml PCR tube
					