



VM **Mini Vortex mixer**



Darsun Scientific Inc, USA
30N Gould St Ste RSheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

FEATURES



Stepless speed regulation
Wide speed range



Brushless motor
Durable



Press to start
Compact body



High quality and
colorful silicone

Darsun Scientific Mini Vortex mixer

DSVM-300/DSVM-300S



Model	DSVM-300	DSVM-300S
Control Mode	Constant speed	Knob
Motor	Brushless motor	
Amplitude	4.5mm (Circular oscillation)	
Maximum Sample Processing Capacity	50ml	
Speed Range	3000rpm	0-3000rpm
Shell Protection Level	IP42	
Power	12W	
Voltage	100-240V, 50/60Hz	
Product Size	93*93*72m	100*93*72mm
Allowable Ambient Temperature	5-40°C, 80%RH	
Gross / Net weight	0.55/0.36kg	0.56/0.37kg