



OS

Overhead Stirrer



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Overhead Stirrer

DSOS-5/DSOS-10/DSOS-15S
/DSOS-30S



FEATURES



Heat more uniformly



304 stainless steel plate,
Corrosion resistant



LCD display



Durable brushless motor

Model	DSOS-5	DSOS-10	DSOS-15S	DSOS-30S
Motor	Brushless motor			
Max Amount Of Stirring	5L	10L	10L	20L
Range Of Rotate Speed	200-1500rpm	100 - 1000rpm	100 - 2000rpm	50 - 1200rpm
Timing	/	/	1min-23h59min	1min-23h59min
Max Torque	10N.cm	20N.cm	150N.cm	30N.cm
Max Viscosity	5000mPas	10000mPas	8000mPas	16000mPas
Drill Chuck Diameter	1-10mm			
Power	20W	25W	30W	36W
Voltage	AC100-240V			