



Premium Vortex Mixer



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Premium Vortex Mixer



Vortex Mixer provides comprehensive mixing capabilities for enzymatic and RIA assays, general test tube mixing, biochemical assays, viral dilution preparation, precipitated assays, cell suspension vortexing and tissue sample mixing.

Specification

- Continuous and Touch modes
- Speed Accuracy: $\pm 5\%$ rpm
- Convenient Speed adjustments by knob control
- Stable for long-term operations
- Multiple accessories for different vessels
- Digital Speed display
- Speed range: 200 to 3000rpm

Model	Premium Vortex Mixer
Speed	Up to 200–3000 rpm
Warranty	2 years
Load Bearing Capacity	0.5 kg
Shipping Weight	14.99 kg
Display	Digital
Dimensions (L x W x H)	290x270x240mm 11.4x10.7x9.5
Motor Type	DC motor
Hertz	50/60 Hz
Material	ABS housing
Includes	Rubber platform
Tube Capacity	≤ 19 tubes
Operating Modes	Touch, Continuous