



COD REACTOR



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

REACTOR

DS-COD-100R



GENERAL FEATURES

- 21 samples at a time.
- Settable heating temperature between 100-165 °C.
- Settable heating time between 0-120 min.
- Transparent plastic protection shield is equipped to protect the operator.

INCLUDES

- 5 pcs turbidity bottles
- Carrying Case

Model	COD-100R	
Measuring Unit	temperature	100-165 °C
	time	0-120 min
General	Power Supply	AC Adapter, 200-240 V AC input
	Dimensions	330 × 228 × 132mm
	Net Weight	4 kg
	Dimensions (Carton)	575X435X285 mm
	Shipping Weight	5.5 kg