



RM **Roller mixer**



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

FEATURES



Timing function



Anti - lock design, more convenient



Adjustable speed



Brushless motor with strong power

Darsun Scientific Roller mixer

DSRM-6H/DSRM-8H



Model	DSRM-6H	DSRM-8H
Motor Type	Brushless motor	
Display	LCD	
Type	Horizontal rolling	
Speed	20-80rpm	
Number Of Rollers	6	8
Roller Length (mm)	280	
Roller Center Distance(mm)	34	
Specs	φ 6-34mm	
Timing	1min-23h59min	
Allowed Ambient Humidity Temp.	5~40℃ 80%RH	
Voltage	AC 100-200V,50?60Hz	
Power	25W	
Dimension (mm)	440*255*105	430*320*100