



Semi-Auto Biochemistry Analyzer,



Darsun Scientific Inc, USA

30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Semi-Auto Biochemistry Analyzer,

DS-100C



Product Features

1. 7"inch true color touch screen(800*480 pixels)
2. 5 USB ports & 1 ethernet port
3. Flash: 256 MB;Over 300 programs, 200000 results can be stored.
4. Resolution: 0.0001 Abs.
5. Photometric range: 0.0000-4.0000 Abs.

Product Overview

Semi-Auto Biochemistry Analyzer is suitable for the analysis of serum, plasma, urine and other body fluid samples, and can complete the test of biochemical, special protein, drug monitoring and blood coagulation items



Multifunctional Incubator With Accurate Constant Temperature
Multifunctional incubator specially designed for chemistry and coagulation tests, different holes can be used to put sample tubes, chemistry cuvettes, coagulation cuvettes, accurate temperature control, more efficient sample incubation



- 5 USB port & Ethernet port
- Mouse, keyboard and printer can be connected. freely easy transmission of is data



Windows embedded compact 7 smart operation system

More convenient to edit test items freely, clear presentation of real-time curve, calibration, quality control and sample results.



7"inch true color touch screen
(800*480) pixels