



HS

Digital Display Hotplate Stirrer



Darsun Scientific Inc, USA
30 N Gould St, Ste R, Sheridan, Wyoming
82801, United States

Email: info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Digital Display Hotplate Stirrer

DSHS-5C/DSHS-5P



FEATURES



Strong magnetic field



Stainless steel ceramic coated countertops

Model	DSHS-5C	DSHS-5P
Control Mode	Brushless motor	Brushless motor
Max Stirring Volume	2L(H ₂ O)	2L(H ₂ O)
RPM	200-1500rpm	50-1500rpm
Timing Range	/	1min-12h
Control Mode	Digital display & Knob	Digital display & Button
Working Plate (mm)	Φ135	Φ135
Plate Temp.	280°C	310°C
Liquid/Plate Temp Accuracy	±1°C/±3°C	±1°C/±3°C
Power	300W	400W
Product Size (mm)	225*155*90	225*155*90