



VM

Vortex Mixer



Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: Info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific



Darsun Scientific

Vortex Mixer

DSM25/DSM25D



- Simple and reliable structure, small dimension, low power consumption and low noise;
- Two models optional: fixed and high/low speed;
- Used for various mixing application with optional adaptor;
- Accept a variety of test sample container: centrifuge tube/PCR plate/cell culture plate/deep well plate;
- Stable operation, small speed error and long service life.

Model	DSM25	DSM25D
Power	30W	30W
Load Adapter	Flat vibrating head	Flat vibrating head / Microplate clamp
Speed Range	2800 rpm	L: 1300 rpm / H: 2800 rpm
Vibration Amplitude	4 mm	
Mix Method	Orbital	
Working Mode	Continuous	
Dimensions (L×W×H)	160 × 160 × 115 mm	
Weight	3.75 kg	
Power Supply	110V/220V; 50/60Hz	