



DS-Mini Transfer (Classical)



Darsun Scientific Inc, USA

30N Gould St Ste R Sheridan, Wyoming
82801, United States

Email: info@darsunscientificinc.com

Phone: +1-3074272810

Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

DS-Mini (Classical) DS-ZYZ10T

Product Features

Simultaneously transfer 1-2 precast or hand-cast gels within an hour. The blue ice box is completely contained within this KIT and absorbs heat generated during quick transfers. The tank of the Transfer Cell can also be used as the tank of our Vertical cell. Perfect match with DS Vertical cell KIT (DS-DYC10T).



Product Overview

The DS Transfer Cell (DS-ZYZ10T) is a part of the Protein Electrophoresis system. It works with the vertical cell (DS- DYC10V) together to provide high-quality blotting of the hand-cast and precast gels. This model accommodates two gel holder cassettes for the transfer of the gels run in the vertical cell (DS-DYC10V).



Model	DS-ZYZ10T
Size	75x100mm
Time	about 60min
Qty of the gels	1-2
Cooling method	Ice Box
Buffer Volume	450ml