



Medical Sanitizer High Pressure Steam Cleaners Machine

Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan,
Wyoming 82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

Medical Sanitizer High Pressure Steam Cleaners Machine

DS – MS8



Product Overview

The Medical Sanitizer High-Pressure Steam Cleaner is designed for **high-temperature steam-based cleaning of medical, surgical, dental, laboratory, and pharmaceutical equipment**. It delivers pressurized steam through a precision hand nozzle to clean surfaces, instrument joints, grooves, trays, and other difficult-to-reach areas.

Product Features

- 1. The use of high temperature and high pressure steam, no need to manually scrub, reduce the labor intensity of staff.
- 2. Suitable for manual pre-cleaning or post-processing of complex and highly polluted instruments, remove stubborn stains attached to the surface of power surgical instruments, ultrasonic knife head and other precision instruments, lumen instruments, ophthalmic instruments, stomatology instruments, obstetrics and gynecology instruments, etc., to achieve the ideal cleaning effect and prevent surface damage of instruments and instruments.
- 3. Continuous steam generation, convenient continuous cleaning, up to 8Bar of high pressure steam can be produced to ensure the cleaning effect of instruments
- .4. Built-in efficient steam generator, designed for precision instrument cleaning in medical environment.

FEATURES



1. Place Instruments



2. Automatic Spray Cleaning



3. Manual Spray Cleaning



4. Final Rinse

DarsunScientific

Technical Specifications

Parameter	Specification
Max. exhaust temperature	160°C
Max. steam pressure	0.8 bar
Pressure chamber volume	6 L
Input Water Type	Automatic/Manual
Voltage	220 V
Power	200 W
Exhaust Air Rate	750 M³/h
Filter liquid collection device	Standard
Voltage	220 V
Power	1200 W

Applications

Hospitals & Healthcare Facilities – Cleaning surgical instruments, trays, racks, and equipment surfaces.

Dental Clinics – Cleaning dental instruments and difficult-to-access areas.

Medical Device Processing – Removing residues from reusable medical equipment before further processing.

Laboratories – Cleaning stainless-steel equipment, work surfaces, trays, and accessories.

Pharmaceutical Industry – Hygienic cleaning of suitable equipment and processing accessories.

Veterinary Clinics – Cleaning reusable instruments and equipment.

Central Sterile Services Departments (CSSD) – Supporting the cleaning stage before sterilization.

Typical Cleaning Process

Instrument Loading → Pre-Rinse → High-Pressure Steam Cleaning → Residue Removal → Rinse/Dry → Inspection → Sterilization