



ULTRASONIC CLEANER

Darsun Scientific Inc, USA
30N Gould St Ste R Sheridan,
Wyoming 82801, United States

Email: Info@darsunscientificinc.com
Phone: +1-3074272810
Web: www.darsunscientificinc.com

DarsunScientific

Darsun Scientific

ULTRASONIC
CLEANER

DS- 150J/225J/300J



Product Overview

Ultrasonic cleaner is a non-polluting, non-damaging, multi-purpose cleaning equipment. It not only has a powerful cleaning effect on objects with complex shapes, but also accelerates dissolution and emulsification. Ultrasonic cleaner has the unparalleled advantages of other physical cleaning or chemical cleaning, and is widely used in service industry, electronics industry, pharmaceutical industry, laboratory, machinery industry, cemented carbide industry, chemical industry and other fields.

Product Feature

- 1. Standard accessories: 1pc of soundproof cover, 1pc of basket, 1pc of PVC pipe.
- 2. The shell and inner cavity are made of 304 stainless steel, which is anti-corrosion.
- 3. LED display for UC-D series.

Technical Parameters

Model	DS-150J	DS-225J	DS-300J
Capacity L	15	22.5	30
Ultrasonic Frequency	Single frequency, 40 kHz	Single frequency, 40 kHz	Single frequency, 40 kHz
Time Range / Heating Temp	1–60 min / RT–80°C	1–60 min / RT–80°C	1–60 min / RT–80°C
Ultrasonic Power (W)	360	480	600
Heating Power (W)	500	500	500
Power Supply	AC 220 V, 50 Hz	AC 220 V, 50 Hz	AC 220 V, 50 Hz
Tank Size (mm)	330 × 300 × 150	500 × 300 × 150	500 × 300 × 200
External Size (mm)	360 × 330 × 310	530 × 330 × 310	530 × 330 × 350
Packing Size (mm)	460 × 400 × 400	640 × 400 × 400	640 × 400 × 400
Gross Weight (kg)	10	13	15